




Milex Pessaries

Pessaries are made of silicone unless otherwise noted and are fully autoclavable.

Pessary Order Information:

- Order by sizes in the medium blue shaded column



- The light blue shade areas show the assortment necessary to fit 85% of your patients.
- Fitting sets available.

To find out more about Pessaries, or to place an order, call 800.243.2974.

Pessary Reimbursement Information:

CPT Codes (silicone)	
Procedure	57160
Supply	99070
HCPCS	A4562

CPT Codes (latex)	
Procedure	57160
Supply	99070
HCPCS	A4561



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611

©2011 CooperSurgical, Inc.

Form #81360 • 04/11

MILEX™ Pessary Product Reference Chart



"We recently used the CooperSurgical Gehrung pessary on a number of patients who underwent hysterectomies years ago and suffer from vault prolapse. These patients prefer nonsurgical treatment and to our delight, the Gehrung works like a charm and my patients are happy and grateful!"

Jonathan Y. Song, MD
TLC Medical Group, St. Charles IL

Ring with Support



KPRS

Size	in.	mm
0	1 3/4"	44
1	2"	51
2	2 1/4"	57
3	2 1/2"	64
4	2 3/4"	70
5	3"	76
6	3 1/4"	83
7	3 1/2"	89
8	3 3/4"	95
9	4"	102
10	4 1/4"	108
11	4 1/2"	114
12	4 3/4"	121
13	5"	127

Shaatz



KPSH

Size in	mm
1 1/2"	38
1 3/4"	44
2"	51
2 1/4"	57
2 1/2"	64
2 3/4"	70
3"	76
3 1/4"	83
3 1/2"	89
3 3/4"	95

Incontinence Ring*



KPCON

Size	in.	mm
0	1 3/4"	44
1	2"	51
2	2 1/4"	57
3	2 1/2"	64
4	2 3/4"	70
5	3"	76
6	3 1/4"	83
7	3 1/2"	89
8	3 3/4"	95
9	4"	102
10	4 1/4"	108
11	4 1/2"	114
12	4 3/4"	121
13	5"	127


Incontinence Dish*



KPCOND

Size mm
55
60
65
70
75
80
85


Gehring*



KPGS

Size	in.(A)	in.(B)	in.(C)	mm(A)	mm(B)	mm(C)
0	1 1/16"	1 1/2"	1 3/8"	27	38	35
1	1 1/8"	1 3/4"	1 1/2"	29	44	38
2	1 3/8"	2"	1 5/8"	35	51	41
3	1 1/2"	2 1/8"	1 3/4"	38	54	44
4	1 5/8"	2 1/4"	1 7/8"	41	57	48
5	1 3/4"	2 3/8"	2"	44	60	51
6	1 7/8"	2 1/2"	2 1/8"	48	64	54
7	2"	2 7/8"	2 1/4"	51	73	57
8	2 1/8"	3"	2 3/8"	54	76	60
9	2 1/4"	3 1/8"	2 1/2"	57	79	64
10	2 3/8"	3 1/4"	2 5/8"	60	83	67

Hodge*



KPEH

Size	in.(A)	in.(B)	mm(A)	mm(B)
0	2 3/4"	1 3/4"	70	44
1	3"	1 7/8"	76	48
2	3 1/4"	1 15/16"	83	49
3	3 1/2"	2"	89	51
4	3 5/8"	2 1/8"	92	54
5	3 3/4"	2 1/4"	95	57
6	3 7/8"	2 3/8"	98	60
7	4 1/4"	2 1/2"	108	64
8	4 5/8"	2 5/8"	117	67
9	5"	2 3/4"	127	70

* Remember for sizing: The incontinence knob available for this pessary adds 1/2 inch (13 mm) to the diameter.

Donut



KPDO

Size in	mm
2"	51
2 1/4"	57
2 1/2"	64
2 3/4"	70
3"	76
3 1/4"	83
3 1/2"	89
3 3/4"	95


Flexible Gellhorn



KPGE

Size in	mm
1 1/2"	38
1 3/4"	44
2"	51
2 1/4"	57
2 1/2"	64
2 3/4"	70
3"	76
3 1/4"	83
3 1/2"	89
3 3/4"	95

Cube**



KPEC

Size	in.	mm
0	1"	25
1	1 3/16"	30
2	1 3/8"	35
3	1 1/2"	38
4	1 5/8"	41
5	1 3/4"	44
6	2"	51
7	2 1/4"	57

** Removed nightly by patient


Inflatoball** (Latex)



KPINF

Size	in.	mm
S	2"	51
M	2 1/4"	57
L	2 1/2"	64
XL	2 3/4"	70


Regula



KPREG

Size	in.(A)	in.(B)	in.(C)	mm(A)	mm(B)	mm(C)
0	1"	1 15/16"	1 3/8"	25	49	35
1	1 1/16"	2"	1 1/2"	27	51	38
2	1 1/8"	2 1/8"	1 5/8"	29	54	41
3	1 1/4"	2 3/16"	1 3/4"	32	56	44
4	1 3/8"	2 1/4"	1 7/8"	35	57	48
5	1 7/16"	2 3/8"	2"	37	60	51
6	1 1/2"	2 5/8"	2 1/8"	38	67	54
7	1 5/8"	2 7/8"	2 1/4"	41	73	57
8	1 3/4"	3 1/8"	2 3/8"	44	79	60

Smith



KPES

Size	in.(A)	in.(B)	mm(A)	mm(B)
0	3 1/8"	2"	79	51
1	3 1/4"	2 1/8"	83	54
2	3 1/2"	2 1/4"	89	57
3	3 3/4"	2 3/8"	95	60
4	4 1/4"	2 1/2"	108	64
5	4 1/2"	2 5/8"	114	67
6	4 3/4"	2 3/4"	121	70
7	5"	2 7/8"	127	73
8	5 1/2"	3"	140	76
9	5 3/4"	3 1/8"	146	79

Order by size as shown in shaded areas.

